

**CPK697Ge11 100µg**  
**BSA Conjugated 3,5-Dinitrosalicylic Acid (DNS)**  
*Instruction manual*

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---


---

11th Edition (Revised in May, 2016)

## [ PROPERTIES ]

**Source:** Protein Conjugation

**Original Structure:**



**Original Chemical Formula:** C<sub>7</sub>H<sub>4</sub>N<sub>2</sub>O<sub>7</sub>

**Original Mol. Mass:** 228.1Da

**Purity:** >90%

**Traits:** Freeze-dried powder.

**Buffer Formulation:** PBS, pH7.4.

**Applications:** Immunogen; Coating Antigen; ELISA; SDS-PAGE.

## [ USAGE ]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## [ STORAGE AND STABILITY ]


**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

## [ IDENTIFICATION ]



**Figure 1. UV-spectrum**