

PAF561Hu61

Cy3-linked Antibody to Uroplakin 2 (UPK2)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PRODUCT INFORMATION]

Immunogen: UPK2, Human

Clonality: Polyclonal

Conjugation: Cy3

Host: Rabbit

Immunoglobulin Type: IgG

Purification: Affinity Chromatography.

Applications: WB, ICC, IHC-P, IHC-F, ELISA

Concentration: 1mg/mL

UOM: 500µg

[IMMUNOGEN INFORMATION]

Immunogen: Recombinant UPK2 (Asp26~Thr153) with two N-terminal Tags, His-tag and T7-tag expressed in *E.coli*.

Accession No.: RPF561Hu01

[ANTIBODY SPECIFICITY]

The antibody is a rabbit polyclonal antibody raised against UPK2. It has been selected for its ability to recognize UPK2 in immunohistochemical staining and western blotting.

[APPLICATIONS]

Western blotting: 1:50-400

Immunocytochemistry in formalin fixed cells: 1:50-500

Immunohistochemistry in formalin fixed frozen section: 1:50-500

Immunohistochemistry in paraffin section: 1:10-100

Enzyme-linked Immunosorbent Assay: 1:100-200

Optimal working dilutions must be determined by end user.


[CONTENTS]

Form & Buffer: Supplied as solution form in PBS, pH7.4, containing 0.02% NaN₃.

[STORAGE]

Store at 4°C for frequent use.

[IMAGES]



Used in Western Blot, Sample:
Recombinant UPK2, Human