

EPB540Ra61 100µg

Eukaryotic Cluster Of Differentiation 147 (CD147)

Organism Species: *Rattus norvegicus* (Rat)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Eukaryotic expression

Host: 293F cell

Residues: Ala22~Arg324

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.9

Predicted Molecular Mass: 34.6kDa

Accurate Molecular Mass: 44-55kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affect the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

```
                AAGFLKAPM  SEQEWAGGSV  VLHCEAVGSP
MPEIQWVFEG  NEPNDSCSQL  WDGARLDRVH  IHATYRQHAA  STLSVDGLAA
EDTGTYECRA  SSSPDRNHLT  RPPRVKVVRA  QASVVVLEPG  TIVTSVQEVD
SKTQLTCFLN  SSGIDIVGHR  WMRGGKVLQE  DTLPLQMKY  TVDADDRSGE
YSCIFLPEPV  GRGNINVEGP  PRIKVGKKSE  HASEGEFVKL  ICKSEASHPP
VDEWVWFKTS  DTGDQTIISNG  TEANSKYVII  STELSELII  SLDLMNVDPG
TYVCNATNSQ  GSARETISLR  VRSR
```

[**IDENTIFICATION**]



Figure. SDS-PAGE

[**IMPORTANT NOTE**]

The kit is designed for research use only, we will not be responsible for any issue if the kit was

used in clinical diagnostic or any other procedures.