



MSI165Cp11

## Medium for Primary Caprine Peripheral Blood Mononuclear Cells (PBMC)

*Instruction manual*

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Revised in Apr, 2025)

### [ Description ]

Medium for Primary Caprine Peripheral Blood Mononuclear Cells, is a liquid medium free of mycoplasma, bacteria and fungi and contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and fetal bovine serum (10%). It is a complete medium designed for optimal growth of Primary Caprine Peripheral Blood Mononuclear Cells *in vitro*.

### [ Components ]

500 ml of Medium for Primary Caprine Peripheral Blood Mononuclear Cells consists of basal medium, 50 ml of FBS, 5 ml of Mononuclear Cell growth supplement and 5 ml of penicillin/streptomycin solution.

*Note: FBS, Mononuclear Cell growth supplement and P/S solution are not pre-mixed in Primary Caprine Peripheral Blood Mononuclear Cells Medium; they must be added separately to make the complete Primary Caprine Peripheral Blood Mononuclear Cells Medium.*

### [ Storage ]

Store the basal medium at 4°C, FBS, Mononuclear Cell growth supplement and P/S solutions at -20°C.

Protect from light.

### [ Shipping ]

Basal medium : room temperature. FBS, Mononuclear Cell growth supplement and P/S solutions: dry ice.

### [ Usage ]

FBS, Mononuclear Cell growth supplement and P/S solutions were thawed at room temperature, sprayed with 75% ethanol on bottles and tubes, and transferred to a sterile operating table. In a sterile field, FBS, Mononuclear Cell growth supplement and P/S solutions were added to the basal medium and mixed evenly to obtain the recombinant complete medium, which was directly used for culture of Primary Caprine Peripheral Blood Mononuclear Cells *in vitro* .

*Note: When stored in the dark at 4°C, the reconstituted complete medium is stable for three months.*

### [ Important note ]

- In order to maintain the best use effect of this product, do not place it in room temperature or high temperature environment for a long time.
- This product is for scientific research use only. It is not for diagnostic, therapeutic, clinical, family and other purposes.