



MSI244Mu11

## Medium for Primary Mouse Placental Mesenchymal Stem Cells (PMSC)

*Instruction manual*

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Revised in Feb, 2024)

### [ Description ]

Medium for Mouse Placental Marrow-derived Mesenchymal Stem Cells, is a liquid medium free of mycoplasma, bacteria and fungi and contains essential and non-essential amino acids, vitamins, organic and inorganic compounds, hormones, growth factors, trace minerals and 10% concentrations of fetal bovine serum . It is a complete medium designed for optimal growth of Mouse Placental Marrow-derived Mesenchymal Stem Cells *in vitro*.

### [ Components ]

500 ml of Medium for Primary Mouse Placental Marrow-derived Mesenchymal Stem Cells consists of basal medium, 50 ml of FBS, 5 ml of Mesenchymal Stem Cell Growth Supplement and 5 ml of penicillin/streptomycin solution.

*Note: FBS, Mesenchymal Stem Cell Growth Supplement and P/S are not pre-mixed in Primary Mouse Placental Marrow-derived Mesenchymal Stem Cells Medium; they must be added separately to make the complete Primary Mouse Placental Marrow-derived Mesenchymal Stem Cells Medium.*

### [ Storage ]

Store the basal medium at 4°C, Mesenchymal Stem Cell Growth Supplement, FBS and P/S solutions at -20°C.

Protect from light.

### [ Shipping ]

Basal medium : room temperature. Mesenchymal Stem Cell Growth Supplement, FBS and P/S solution: dry ice.

### [ Usage ]

Mesenchymal Stem Cell Growth Supplement, FBS and P/S solutions were thawed at room temperature, sprayed with 75% ethanol on bottles and tubes, and transferred to a sterile operating table. In a sterile field, Mesenchymal Stem Cell Growth Supplement, FBS and P/S solutions were added to the basal medium and mixed evenly to obtain the recombinant complete medium, which was directly used for culture of Primary Mouse Placental Marrow-derived Mesenchymal Stem Cells *in vitro* .

*Note: When stored in the dark at 4°C, the reconstituted complete medium is stable for three months.*

### [ Important note ]

- In order to maintain the best use effect of this product, do not place it in room temperature or high temperature environment for a long time.



- This product is for scientific research use only. It is not for diagnostic, therapeutic, clinical, family and other purposes.