

CPA601Bo21 200µg
OVA Conjugated Myeloperoxidase (MPO)
Organism Species: *Bos taurus*; Bovine (Cattle)
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Protein Conjugation

Residues: Thr38-Glu53

Purity: >90%

Traits: Freeze-dried powder

Buffer Formulation: 10mM PBS, pH7.4.

Predicted isoelectric point: 5.9

Predicted Molecular Mass: 1942.3Da

Applications: Immunogen; Blocking Peptide.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

TCM EEA~~K~~RLVDKV YKE

[IDENTIFICATION]




Figure 1. MASS Spectrometry of the Original Peptide




Figure 2. SDS-PAGE of OVA Conjugated MPO

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.