



CSI282Ra01

Primary Rat Dorsal Root Ganglion Neuron Cells (DRGN)

Organism Species: *Rattus norvegicus* (Rat)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

1st Edition (Revised in May, 2024)

[DESCRIPTION]

Cell Type: Neuron Cell

Synonyms: DRGN

Species: *Rattus norvegicus* (Rat)

Tissue Source: Dorsal root ganglion

Size: $>5 \times 10^5$ cell/vial

[PROPERTIES]

Cell activity: $>85\%$ (Viability by Trypan Blue Exclusion).

Formulation: Frozen 1 mL or T25 flask.

Biosafety: Negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast and fungi.

Applications: For research use only. It is not approved for human or animal use, or for application in clinical diagnostic procedures.

Growth Properties: Adherent

[CONTENTS]

Form & Buffer: Supplied as solution form in frozen stock solution, containing 90% FBS+10% DMSO.

[USAGE]

Upon receiving the cells in a T-25 flask at room temperature, immediately transfer the cells to 37°C, 5% incubator; the cells in vials, directly and immediately transfer the cells from dry ice to liquid nitrogen.

Culture conditions:

Special culture medium for neuronal cell:

Neurobasal-A Medium+B-27 Supplement (50X)+1%Penicillin-Streptomycin Solution

Temperature: 37°C

Condition: 95% air, 5% carbon dioxide

Cell recovery:

After receiving the cells, shake at 37°C in a water bath until completely dissolved, transfer to a 15 ml centrifuge tube, add 3-5 times complete culture solution, 1000 rpm for 5 min, discard the supernatant, and place in a T25 flask for culture.

Cell passage:

Further culture of Rat Brain Stem Neuron Cell is guaranteed under the conditions we provide; however,



it is not recommended for expansion or long-term culture because is a terminally differentiated cell and does not proliferate in culture.

[Shipping]

Dry ice.

[STORAGE]

Upon receiving, directly and immediately transfer the cells from dry ice to liquid nitrogen and keep the cells in liquid nitrogen until they are needed for experiments.

[IMPORTANTNOTE]

1. Rat Dorsal Root Ganglion Neuron Cells is not recommended for expanding or long-term cultures since the cells do not proliferate in culture.
2. The cell is for research use only, we will not be responsible for any issue if the cell was used in clinical diagnostic or any other procedures.

[Figure]




Figure 1

Figure 1 Morphology of Rat Dorsal Root Ganglion Neuron Cells (Optical microscope, ×100)




Figure 2

Figure 2 Morphology of Rat Dorsal Root Ganglion Neuron Cells (Optical microscope, ×200)




Figure 3

Figure 3 Immunofluorescence identification of Tubulin Beta (TUBb) specific antibody (×200)




Figure 4

Figure 4 Immunofluorescence identification of Tubulin Beta (TUBb) specific antibody (×400)