

**RPD686Ra01 50µg**  
**Recombinant 3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase (HMGCR)**  
**Organism Species: Rattus norvegicus (Rat)**  
***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

---

12th Edition (Revised in Aug, 2016)

## **[ PROPERTIES ]**

**Source:** Prokaryotic expression.

**Host:** *E. coli*

**Residues:** Ala638~Ala887

**Tags:** N-terminal His-Tag

**Subcellular Location:** Endoplasmic reticulum membrane; Multi-pass membrane protein.

**Purity:** >99%

**Traits:** Freeze-dried powder

**Buffer formulation:** 100mM NaHCO<sub>3</sub>, 500mM NaCl, pH8.3, containing 1mM EDTA, 1mM DTT, 0.01% sarcosyl, 5%Trehalose and Proclin300.

**Original Concentration:** 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 8.4

**Predicted Molecular Mass:** 30.0kDa

**Accurate Molecular Mass:** 30kDa as determined by SDS-PAGE reducing conditions.

## **[ USAGE ]**

Reconstitute in 100mM NaHCO<sub>3</sub>, 500mM NaCl (pH8.3) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## **[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ SEQUENCE ]**

```
                                AGR NLYIRLQSKT
GDAMGMNMIS KGTEKALLKL QEFFPELQIL AVSGNYCTDK KPAAINWIEG
RGKTVVCEAV IPAKVVREVL KTTTEAMVDV NINKNLVGSA MAGSIGGYNA
HAANIVTAIY IACGQDAAQN VGSSNCITLM EASGPTNEDL YISCTMPSIE
IGTVGGGTNL LPQQACLQML GVQGACKDNP GENARQLARI VCGTVMAGEL
SLMAALAAGH LVRSHMVHNR SKINLQDLQG TCTKKA
```

**[ IDENTIFICATION ]**




Figure 1. SDS-PAGE