

RPN267Hu01 50µg

Recombinant Ankyrin Repeat And KH Domain Containing Protein 1 (ANKHD1)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Val204~Ala563

Tags: N-terminal His Tag

Subcellular Location: Cytoplasm

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 50µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 4.2

Predicted Molecular Mass: 41.8kDa

Accurate Molecular Mass: 43kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

VDTRSLA EACSDGDVNA VRKLLDEGRS VNEHTEEGES LLCLACSAGY
YELAQVLLAM HANVEDRGNK GDITPLMAAS SGGYLDIVKL LLLHDADVNS
QSATGNTALT YACAGGFVDI VKVLLNEGAN IEDHNENGHT PLMEAASAGH
VEVARVLLDH GAGINTHSNE FKESALTLAC YKGHLDMVRF LLEAGADQEH
KTDEMHTALM EACMDGHVEV ARLLLDGSAQ VNMPADSFES PLTLAACGGH
VELAALLIER GANLEEVNDE GYTPLMEAR EGHEEMVALL LAQGANINAQ
TEETQETALT LACCGGFSEV ADFLIKAGAD IELGCSTPLM EASQEGHLEL
VKYLLASGAN VHA
    
```

[IDENTIFICATION]

GGATCGGAGCCTCCAGCTCGGAGCCTCTGTGAGTGGGGTGTATACCTCTGTGAAATCTTAGTGGAGGAGAGTGTGAGATGAGGAGGGGAGGAGGCTCTGTGTTTCTCTTCCAGAGGATTAAGATTCAGAGATTTCTTCTTAGGATGCTATGTGTGAGATGGAGGATTAACGAGATTAATCTCCCGAGCCCTCCAGCGGCTTCTAGATTTGGAAAT




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.