

RPD546Mu01 100ug

Recombinant Cytochrome P450 11B1 (CYP11B1)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Asn501

Tags: Two N-terminal Tags, His-tag and SUMO-tag

Subcellular Location: Mitochondrion

Purity: > 80%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 9.8

Predicted Molecular Mass: 71.4kDa

Accurate Molecular Mass: 71kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.


Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MTMALRVTTDVWLARPWQCLHRTRALGTTATLAPKTLQPFEAIPQYSRNKWLKMIQILREQ
GQENLHLEMHQVFRELGPFRHSVKGKTQIVFVTLPEDVEKLYQVESTHPCRMPLESWIVHRE
LRGLGRGVFLLNGPEWYFNRLQLNPNVLSKAVQKQFVPLVDGIARDFVDNLKKKMLESVHG
SFSMDFQSSVFNYTIEASHFVLFGERLGLIGRDLSPDSLKFLHTLHSMFKTTTQLLYLPRSLTR
WTSTRVWKENLESWDFISEYVTKCIKNVYRELAEGRPQSWSVTAELVAERTLSMDAIQANS
MELIAGSTDTTSTPLVMTFFELARNPDVQQALRQESLAAEASIAANPQKAMSDLPLLRALK
ETRLYPVGTFLERILSSDLVLQNYHVPAGTVLNVNLYSMGRNPAVFPRPERYMPQRWLERK
RSFKHLAFGFGVRQCLGRRLAEAEEMMLLLHHVLSFHVETQEKEDVRMAYRFVLMPPSSSPL
LTFRPVN

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.