

RPB240Mu01 100µg

Recombinant Cytokeratin 20 (CK 20)

Organism Species: *Mus musculus* (Mouse)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Asp2~Val431

Tags: N-terminal His Tag

Subcellular Location: Cytoplasm

Purity: > 97%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.0

Predicted Molecular Mass: 50.2kDa

Accurate Molecular Mass: 50kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

DFSRQSFHR SLSSSSQGPA LSMGSLYRK GTVQRLGAAP SVYGGAGGHG
TRISVSKAVM SYGGDLSNGS DLFGGNGKLA MQNLNDRLAN YLEKVRSEQ
SNSRLEAIK QWYETNAPST IRDYSSYYAQ IKELQNQVKD AOVQNAQCVL
RIDNAKLAAE DFRLKFETER GMRIAVEADL QGLSKVYDNL TLQKTDLEIQ
IEELNKDLAL LKKEHQEEVE VLRRQLGNNV NVEVDAAPGL NLGEIMNEMR
QRYEVLAKN LQEAKEQFER QSOTLQQQVT VNTEELKGFE VQVTELRRTY
QNLEIELQSH LSMKESLERN LEDVKARYAS QLAAIQEMLS SLEAQLMQIR
SDTERQNQEH NILLDIKTRL EQEIATYRRL LEGEDIKTTE YQLSTLEMKD
IKKTRKIKTV VEEVVDGKVV SSEVKEIEES V
    
```

[IDENTIFICATION]




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.