



**RPB310Mu01 10 $\mu$ g**  
**Recombinant Elastase 4 (ELA4)**  
**Organism Species: *Mus musculus* (Mouse)**  
***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

---

12th Edition (Revised in Aug, 2016)

## [ **PROPERTIES** ]

**Source:** Prokaryotic expression.

**Host:** *E. coli*

**Residues:** Val30~Leu268

**Tags:** N-terminal His-Tag

**Tissue Specificity:** Pancreas.

**Purity:** >98%

**Traits:** Freeze-dried powder

**Buffer formulation:** 10mM PBS, pH7.4, containing 1mM DTT, 5% trehalose, 0.01% sarcosyl and Proclin300.

**Original Concentration:** 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 5.8

**Predicted Molecular Mass:** 30.3kDa

**Accurate Molecular Mass:** 30kDa as determined by SDS-PAGE reducing conditions.

## [ **USAGE** ]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## [ **STORAGE AND STABILITY** ]

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.


**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ SEQUENCE ]**

```

                                V VGGEDAVPNS WPWQVSLQYL
RDDTWRHTCG GSLITTSHVL TAAHCINTNL TYRVGLGKYN LTVEDEEGSV
YAEVDTIYVH EKWNRLLLWN DIAIIKLAEP VELSDTIQVA CIPEQDSLLP
GDYPCYVTGW GRLWTNGPIA EVLQQGLQPI VNHTTCSRLD WWFIKVRETM
VCAGGDGVIS ACNGDSGGPL NCPVEDGLWQ VHGIVSFGSS RGCNTYKKPV
VFTRVSAYID WIKEKIQL
```

**[ IDENTIFICATION ]**



**Figure 1. SDS-PAGE**