

RPA337Hu01 100µg

Recombinant GATA Binding Protein 2 (GATA2)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Pro175~Gly480

Tags: N-terminal His Tag

Subcellular Location: Nucleus

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose .

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 9.7

Predicted Molecular Mass: 36.2kDa

Accurate Molecular Mass: 36kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.


Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

VSPDPSTTGA	ASPASSSAGG	SAARGEDKDG	VKYQVSLTES	MKMESGSPLR	PGLATMGTQP	PTPPKE
ATHHPIPTYF	SYVPAAAHDY	SSGLFHPGGF	LGGPASSFTP	KQRSKARSCS	EGRECVNCGA	
TATPLWRRDG	TGHYLCNACG	LYHKMNGQNR	PLIKPKRRLS	AARRAGTCCA	NCQTTTTTLW	
RRNANGDPVC	NACGLYYKLH	NVNRPLTMKK	EGIQTRNRKM	SNKSKKSKKG	AECFEELSKC	
MQEKSSPFA	AALAGHMAPV	GHLPPFSHSG	HILPTPTPIH	PSSSLFSGHP	HPSSMVTAMG	

[IDENTIFICATION]



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.