

RPE908Hu01 10µg
Recombinant Histone Deacetylase 8 (HDAC8)
Organism Species: *Homo sapiens (Human)*
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Val377

Tags: N-terminal His Tag

Subcellular Location: Nucleus

Purity: > 80%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 100µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.4

Predicted Molecular Mass: 45.5kDa

Accurate Molecular Mass: 46&22kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MEEPEEPADS GQSLVPVYIY SPEYVSMCDS LAKIPKRASM VHSLIEAYAL HKQMRIVKP
KVASMEEMAT FHTDAYLQHL QKVSQEGDDD HPDSIEYGLG YDCPATEGI FDYAAAIGGA
TITAAQCLID GMCKVAINWS GGWHHAKKDE ASGFCYLND AVLGILRLRR KFERILYVDL
DLHHGDGVED AFSFTSKVMT VSLHKFSPG FFPGTGDVSD VGLGKGGRYS VNVPIQDGIQ
DEKYYQICES VLKEVYQAF NPKAVVLQLG ADTIAGDPMC SFNMTPVGIG KCLKYILQWQ
LATLILGGG GYNLANTARC WTYLTGVILG KTLSSEIPDH EFFTAYGPDY VLEITPSCR
PDRNEPHRIQ QILNYIKGNL KHVV

[IDENTIFICATION]




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.