

RPB983Ra01 200µg

Recombinant Microtubule Associated Protein Tau (MAPt)

Organism Species: *Rattus norvegicus* (Rat)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Ala34~Ser368

Tags: N-terminal His Tag

Subcellular Location: Membrane, Cytoplasm

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 9.9

Predicted Molecular Mass: 38.4kDa

Accurate Molecular Mass: 48kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affect the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

```
                                AEEAGIG DTPNMEDQAA
GHVTQARVAG VSKDRTGND E KKAKGADGKT GAKIATPRGA ATPGQKGTSN
ATRIPAKTTP SPKTPPGSGE PPKSGERSGY SSPGSPGTPG SRSRTPSLPT
PPTREPKKVA VVRTPPKSPS ASKSRLQTAP VPMPDLKNVR SKIGSTENLK
HQPGGGKVQI INKKLDLSNV QSKCGSKDNI KHVPGGGSVH IVYKPVDLSK
VTSKCGSLGN IHKPGGGQV EVKSEKLD FK DRVQSKIGSL DNITHVPGGG
NKKIETHKLT FRENAKAKTD HGAEIVYKSP VVSGDTSPRH LSNVSSTGSI
DMVDSPQLAT LADEV SAS
```

[**IDENTIFICATION**]

```
JGCTGAGAGACAGCCATCGGAGACACCCCGACATCGAGGACCAAGCTCTGGCCATGTGACTCAAGCTGAGTGGCCGCGTAAAGCAAGGACAGGACAGSAAATGACGAGAGAGACCCAGGCGGCGGATGCGAAGAGCGGCGGAGGATCGCCACCTCGGAGGCAAGCCACTCGGGGCGAGAAAGGCGCATCCATTCGACCGA
A E E A G I G D T P N H E D Q A A G H V T Q A R V A G V S K D R T G N D E K K A K G A D G K T G A K I A T P R G A A T P G Q K G T S N A T R
```




Figure. Gene Sequencing (Extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.