

RPC250Hu01 100µg

Recombinant Pirin (PIR)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[**PROPERTIES**]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Met1~Asn290

Tags: N-terminal His Tag

Subcellular Location: Nucleus, Cytoplasm

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: 100mMNaHCO₃, 500mMNaCl, pH8.3, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 500µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.4

Predicted Molecular Mass: 35.8KDa

Accurate Molecular Mass: 36kDa as determined by SDS-PAGE reducing conditions.

[**USAGE**]

Reconstitute in ddH₂O to a concentration of 0.1-0.5 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

MGSSKKVTL S VLSREQSEGV GARVRRS IGR PELKNLDPFL LFDEFKGRP
 GGFDPHPRG FETVSYLLEG GSM AHEDFCG HTGKMNPGDL QWMTAGRGIL
 HAEMPCSEEP AHGLQLWVNL RSSEKMVEPQ YQELKSEEIP KPSKDGVTVA
 VISGEALGIK SKVYTRPTL YLDFKLDPGA KHSQPIPKGW TSFIYTISGD
 VYIGPDDAQQ KIEPHHTAVL GEGDSVQVEN KDKRSHFVL IAGEPLREPV
 IQHGPFVMNT NEEISQAILD FRNAKNGFER AKTWKSKIGN

[IDENTIFICATION]

GAATTGATGGGCTCTCCAGAAAGTTACTCTCTGAGTCTGAGCCGGAGCGTGGAGGGGTTGGAGGGGGTGGAGAGCATGGGAGCCGGAGTTAAAGATCTGGATGGGTTTTACTGTTGATGATTTAAAGAGGTAGCCGGAGGATTTCTGATCATCCATCGAGGTTTTGAAAGGATCTCTCTCCAGGGGGC





Figure . Gene Sequencing (extract)



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.