

**RPM277Hu01 100µg**

**Recombinant Testican 2 (TIC2)**

**Organism Species: *Homo sapiens (Human)***

***Instruction manual***

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

13th Edition (Revised in Aug, 2023)

**[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Glu23~Trp424

**Tags:** N-terminal His Tag

**Subcellular Location:** Secreted

**Purity:** > 90%

**Traits:** Freeze-dried powder

**Buffer formulation:** PBS, pH7.4, containing 0.01% SKL, 5% Trehalose .

**Original Concentration:** 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 4.4

**Predicted Molecular Mass:** 48.3kDa

**Accurate Molecular Mass:** 48kDa as determined by SDS-PAGE reducing conditions.

**[ USAGE ]**

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

**[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ SEQUENCE ]**

EGDAKGLK EGETPGNFME DEQWLSSISQ  
YSGKIKHWNR FRDEVEDDYI KSWEDNQGD EALDTTKDPC QKVKCSRHKV  
CIAQGYQRAM CISRKKLEHR IKQPTVKLHG NKDSICKPCH MAQLASVCGS  
DGHTYSSVCK LEQQACLSSK QLAVRCEGPC PCPTEQAATS TADGKPETCT  
GQDLADLGDR LRDWFQLLHE NSKQNGSASS VAGPASGLDK SLGASCKDSI  
GWMFSKLDTS ADLFLDQTEL AAINLDKYEV CIRFFFNSCD TYKDGRVSTA  
EWCFCFWREK PPCLAELERI QIQEAAKKKP GIFIPSCDED GYYRKMQCDQ  
SSGDCWCVDQ LGLELTGTRT HGSPDCDDIV GFSGDFGSGV GWEDEEEKET  
EEAGEEAEFE EGEAGEADDG GYIW

**[ IDENTIFICATION ]**




Figure. SDS-PAGE

**[ IMPORTANT NOTE ]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.