

RPA950Hu01 1mg
Recombinant Von Willebrand Factor Cleaving Protease (vWFCP)
Organism Species: *Homo sapiens (Human)*
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Val896~Val1131

Tags: N-terminal His and GST Tag

Subcellular Location: Secreted

Purity: > 80%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose .

Original Concentration: 150µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.6

Predicted Molecular Mass: 55.4kDa

Accurate Molecular Mass: 60kDa as determined by SDS-PAGE reducing conditions.

Phenomenon explanation:

The possible reasons that the actual band size differs from the predicted are as follows:

1. Splice variants: Alternative splicing may create different sized proteins from the same gene.
2. Relative charge: The composition of amino acids may affects the charge of the protein.
3. Post-translational modification: Phosphorylation, glycosylation, methylation etc.
4. Post-translation cleavage: Many proteins are synthesized as pro-proteins, and then cleaved to give the active form.
5. Polymerization of the target protein: Dimerization, multimerization etc.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

                                VWTPA
AGSCSVSCGR GLMELRFLCM DSALRVPVQE ELCGLASKPG SRREVCQAVP
CPARWQYKLA ACSVSCGRGV VRRILYCARA HGEDDGEEIL LDTQCQGLPR
PEPQEACSL E PCPPRWKVMS LGPCSASCGL GTARRSVACV QLDQGQDVEV
DEAACAALVR PEASVPCLIA DCTYRWHVGT WMECSVSCGD GIQRRRDTCL
GPQAQAPVPA DFCQHLPKPV TVRGCWAGPC V
    
```

[IDENTIFICATION]

G A A T T C G T G T G G A C C C C T G C G G C A G G G T C G T G C T C C G T C T C C T G C G G G C G A G G T C T G A T G G A G C T G C G T T T C C T G T G C A T G G A C T C T G C C C T C A G G G T G C C T G T C C A G G A A G A G C T G T G T G G C C T G G C A




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.